



ProAxis®

Adjust your patient, not your procedure.




Precise Intraoperative Control Never Looked so Good.


Increased Access and Alignment in a L4-L5 TLIF with
Posterior Column Osteotomy (PCO)



ProAxis was flexed to 25°, enabling spine flexion and optimal visualization.




Lordosis increased to 28° following PCO completion and TLIF cage placement.




Gradual extension of ProAxis to 0° enabled controlled closure and additional lordosis.

Preoperative Radiograph



Postoperative Radiograph




L4-L5 lordosis increased 9° from 24° preoperatively to 33° postoperatively.

Experience a new degree

From Simple Degenerative to Complex Deformity Cases

As the first of its kind, ProAxis® is a state-of-the-art, hinged spinal table with multi-procedural versatility. By integrating a 55-degree range of robotically-assisted flexion and extension into a single table, ProAxis enables precise patient positioning and optimal alignment across a wide range of cases.

Dual-mode hinge offers dynamic articulation and precise intraoperative control for unparalleled usability across prone, supine, and lateral procedures.



of versatile access and prec

Targeted Access and Visualization

ProAxis seamlessly delivers up to 35° of flexion to help improve access and visualization at the surgical site.

- **Robotically-assisted hinge** allows for controlled micro-movements and real-time adjustments
- **Intraoperative adjustments** provide enhanced surgical access and direct visualization
- **Radiolucent carbon-fiber construction** enables unobstructed anterior, posterior, lateral and oblique intraoperative imaging
- **Radiolucent open frame and removable two-piece flat-top** provide full compatibility with O-arm® and C-arm imaging systems



Radiolucent hinge enables continuous micro-movement precision.

More Lordosis, Less Stress

With up to 20° of controlled extension, ProAxis enables physicians to optimize lordosis with less force and reduced manual compression on instrumentation and bone.¹



Extending the ProAxis facilitates gradual closure and precise lordotic adjustments.

ise alignment.

Designed with Simplicity in Mind and in Hand

Thoughtfully designed features enhance ease-of-use and workflow while supporting intraoperative patient positioning adjustments.

- **Torso Trolley®** provides translation of the patient's torso during intraoperative flexion and extension to prevent skin shearing
- **Fixed Surgical Site mode** maintains the surgical site position during hinge articulation
- **Handheld IntelliPendant™** provides intuitive controls for table height, hinge angle, lateral tilt, and Trendelenburg position

Torso Trolley features an adjustable Chest Slide and 5-way adjustable arm boards for optimal patient positioning and synchronized movement.



IntelliPendant™ delivers intraoperative control and real-time table status— all at the touch of a button.



The ProAxis Spinal Surgery Table.

Position with precision. Align with accuracy.

The revolutionary ProAxis® Spinal Surgery Table delivers a 55-degree range of robotically-assisted patient positioning to help you achieve optimal alignment, and reach procedural goals without compromise.



Prone

Supine

Lateral

Ease into a new level of precise positioning.

Call **1-800-777-4674** for more information about the ProAxis.



Patient Positioning Area Specifications

Prone	19 in. (48 cm) Width
	81 in. (206 cm) Length

Supine / Lateral	21 in. (53 cm) Width
	81 in. (206 cm) Length

Overall Table Length*	122–138 in. (310–351 cm)
-----------------------	--------------------------

Retracted Storage Length	80 in. (203 cm)
--------------------------	-----------------

Height Range*	20–47 in. (51–119 cm)
---------------	-----------------------

Hinge Up	35°
----------	-----

Hinge Down	20°
------------	-----

Lateral Tilt	Tilt 15° (Right / Left)
--------------	-------------------------


Trendelenburg†	15°
----------------	-----

Reverse Trendelenburg†	20°
------------------------	-----

Patient Weight Capacity	500 lbs. (227 kg)
-------------------------	-------------------

*Varies depending on table's position and motion

†From Home Position



Patent Numbers: US 8,584,281. US 9,958,313.

ProAxis is designed, developed, and manufactured in the USA.

Ordering Information

6988 ProAxis Spinal Surgery Table	100/120 V, 50/60 Hz
-----------------------------------	---------------------

6988I ProAxis Spinal Surgery Table	220/240 V, 50/60 Hz
------------------------------------	---------------------



30031 Ahern Avenue

Union City, CA 94587-1234 USA

Toll Free: 800-777-4674 | Outside USA: +1-510-429-1500

Fax: 510-429-8500

mizuho.com

© 2022 Mizuho OSI. NW1155 Rev A

1. Vedantam A, Verla T, Mayer RR, Raber MR, Ropper AE. Use of an open-frame hinged surgical table to restore segmental lumbar lordosis after posterior column osteotomy. *Int J Spine Surg.* 2020;14(3):316-320.

O-arm® is a registered trademark of Medtronic, Inc.