







Helmadrasser / Mattress Paddings

Text	Item No	Cover Colour	Size in cm LxWxH	
Alle madrasser laves i de herunder tre viste varianter i standardhøjde 8 og 10 cm. Alle er baseret på en højdensitets HR skum i bunden og to typer viskoelastisk skum, der giver en gradueret overgang fra en relativt blød top til en fast bund. Det sikrer at patienten ligger komfortabelt, godt aflastet uden risiko for at ligge igennem madrassen. All mattresses are made in the three variants shown below in standard heights of 8 and 10 cm. All are based on a high-density HR foam in the base and two types of viscoelastic foam that provide a gradual transition from a relatively soft top to a firm base. This ensures that the patient lies comfortably, well-relieved without the risk of lying through the mattress.				
 Pediatric  Standard  Bariatric				
Herunder finder I enkle eksempler på forskellige madrasser Below you will find a few examples of different types of mattresses.				
Helmadras 1-Part mattress padding				
Helmadras 1-Part mattress padding	1111.5041 1114.1206 1111.7709 1111.0175	Black Black	181x53x8 200x51x8 220x52x8 230x54x8	
2-delt madras 2-Part mattress padding				



Helmadrasser / Mattress Paddings

Text	Item No	Cover Colour	Size in cm LxWxH	
2-Delt madras 2-Part mattress padding	1111.7703 1111.7904 1111.7906	Black Black		
Vi har gennem årene lavet et utal forskellige madrasser til operationslejer. I stedet for at opliste alle disse , som måske alligevel ikke indeholder den helt rigtige for jer, har vi valgt at lave en række formularer, som findes på de næste sider, som I kan udfylde og maile til os. Vi vil så matche jeres forespørgsel med allerede eksisterende konfigurationer, finde den der passer eller lave et nyt design og udarbejde et tilbud. Over the years, we have made countless different mattresses for surgical tables. Instead of listing all of these, which may or may not contain the right one for you, we have chosen to create a series of forms on the following pages that you can fill out and email to us. We will then match your enquiry with existing configurations, find the one that fits or create a new design and prepare a quote.				