






TopMadrasser / Overlay Paddings

Text	Item No	Cover Colour	Size in cm LxWxH	
Alle topmadrasser laves i de herunder tre viste varianter i standardhøjde 5 og 7 cm. Alle er baseret på en højdensitets HR skum i bunden og to typer viskoelastisk skum, der giver en gradueret overgang fra en relativt blød top til en fast bund. Det sikrer at patienten ligger komfortabelt, godt aflastet uden risiko for at ligge igennem madrassen. All overlay paddings are made in the three variants shown below in standard heights of 8 and 10 cm. All are based on a high-density HR foam in the base and two types of viscoelastic foam that provide a gradual transition from a relatively soft top to a firm base. This ensures that the patient lies comfortably, well-relieved without the risk of lying through the mattress.				
 Pediatric  Standard  Bariatric				
Herunder finder I enkle eksempler på forskellige topmadrasser Below you will find a few examples of different types of overlay paddings.				
TopMadras 1-Part Overlay padding				
Helmadras 1-Part mattress padding	1120.1005 1113.0119 1114.1001	Black Black	190x51x7 189x52x7 187x54x7	



TopMadrasser / Overlay Paddings

Text	Item No	Cover Colour	Size in cm LxWxH	
Vi har gennem årene lavet et utal forskellige topmadrasser til operationslejer. I stedet for at opliste alle disse , som måske alligevel ikke indeholder den helt rigtige for jer, har vi valgt at lavel en række formularer, som findes på de næste sider, som I kan udfylde og maile til os. Vi vil så matche jeres forespørgsel med allerede eksisterende konfigurationer, finde den der passer eller lave et nyt design og udarbejde et tilbud. Over the years, we have made countless different overlay paddings for surgical tables. Instead of listing all of these, which may or may not contain the right one for you, we have chosen to create a series of forms on the following pages that you can fill out and email to us. We will then match your enquiry with existing configurations, find the one that fits or create a new design and prepare a quote.				