



CONDOR[®] MedTec
EXPAND YOUR POSSIBILITIES



Condor GoldLine[®] for Minimally Invasive Surgery



MITTELSTAND
GLOBAL
HEALTH MADE
IN GERMANY



WWW.CONDOR-MEDTEC.DE/EN



Minimally Invasive Surgery

Laparoscopic Liver Retractor

-  Constant, smooth lifting of the liver
-  Calm situs with optimal view of the upper abdominal organs
-  Flexible application possibilities, such as interventions with robotic system
-  Individual positioning due to ball joints
-  Reusable




Laparoscopic Liver Retractor 100.062.175 and 100.091.212

The ideal extension for a calm situs.

The laparoscopic retractor serves the smooth raising of the liver and allows for an excellent view of the upper abdominal organs.

Placed in a sterile area at the operating table, these retractors may be individually adjusted to the anatomic conditions of the patient thanks to their ball joint.

Available in various sizes. May also be used for obese patients.







Laparoscopic Liver Retractor for Robotics 100.091.430

Ideal for use with a robotic system.

Do you also use one of the robotic arms, to keep the liver out of the field of view?

+ With the GoldLine® liver retractor for robotics you regain the flexible application possibilities of the arm.

At the beginning of surgery, the retractor is placed centrally under the liver and continuously keeps it out of your field of vision in situs. Due to the lateral positioning, both the robotic system and the surgical staff maintain their complete freedom of movement.



Instrument Fixation

For holding instruments or optics during minimally invasive procedures.

The holder for the liver paddle allows you to hold your liver retractors in a stable position on its handle. This gives you maximum length for using the instrument in the situs. Accordingly, the use is ensured even with very obese patients.

- +** Hold your instruments and optics in a stable position.
- +** Adjustable with one hand or firmly fixed thanks to locking hole and ball joint.

100.020.045

Holder for liver paddle
115.110.035

8 mm 115.008.035	12 mm 115.012.035	5 mm 115.005.014	115.014.035	10 mm 115.010.014
---------------------	----------------------	---------------------	-------------	----------------------

Instrument fixation ball

Instrument fixation especially suitable for optics (consisting of base and sleeve).




Laparoscopic Instrument Holder

Benefits and advantages of the support system



Advanced comfort for the surgeon

- Ergonomic working sparing the back
- Store your instruments close to the surgical area
- Prevention of tangled cables
- Optimized process flow by quick change of instruments and a calmer working – developed by users for users



Patient Safety

- No instruments rest upon the patient
- No damaging of the sterile cover by hot or sharp instruments



Compatibility

- Ideal for the use with the Condor® GoldLine liver retractor
- May be extended as a supplement to the existing portfolio (compatible with the standard main holder)



made in
Germany



**Give us a call and arrange
a non-binding trial!**



**Excellent
Service**



**Quality Work
"Made in Germany"**



**Tried
and Tested**